

Birth of Craftsmanship

These original forceps have been created by combining the many years of experience

The forceps are even smaller and more precise for outstanding ease of use.

Through the utmost care and attention, we are proud to put these superior instruments in the hands of users.




ControMed



SUPERMICRO SERIES

Precise and ultra-small. Ultimate form designed

Ear forceps




SuperMicro ear forceps
serrated jaws
Tip: 0.6 (W) × 2.0 (D) mm



Micro ear forceps
serrated jaws, straight
Tip: 1.0 (W) × 4.2 (D) mm

Ear forceps



SuperMicro ear forceps
through-cutting
Internal dimension: 0.6 (W) × 2.8 (D) mm



Takahara ear forceps
through-cutting, straight
Internal dimension: 1.0 (W) × 3.0 (D) mm

Scissors




SuperMicro scissors
Edge: 2.5 mm




Shea micro scissors
straight
Edge: 6.0 mm

Cupped forceps




SuperMicro cupped forceps
oval cup
Tip: 0.6 (W) × 0.8 (D) mm




Micro cupped forceps
straight
Tip: 1.0 (W) × 2.0 (D) mm

Cupped forceps



SuperMicro cupped forceps
round cup
Tip: 0.8 (W) × 0.8 (D) mm

Malleus nipper



SuperMicro
malleus nipper
Edge: 0.8 mm

Takahara forceps



Takahara SuperMicro
cupped forceps
Tip: 0.6 (W) × 2.0 (D) mm



Takahara ear forceps
oval cup, straight
Tip: 1.2 (W) × 4.0 (D) mm

ControMed

POX-BOX: 1056, 26 Walnut Road,
Melbourne, VIC 3104, Australia
Tel : +61 4 3074 2079
Fax: +61 3 9891 7050
Website: www.contromed.com